



Your Vision, Our Future

Date: 07/17/2015

The following list of Microsoft **Windows XP** updates have been tested and approved for **EndoWorks 7.4.42** compatibility. Prior to applying Windows Updates, make sure your system is current with EndoWorks patches. For questions regarding applying Microsoft Updates or EndoWorks Patches, please contact Olympus Technical Assistance Center.

NOTE: Do not update Internet Explorer (IE) to a version other than the one supported on your EndoWorks system. Refer to the EndoWorks 7.4 IT Specifications (TR0020) document for a list IE versions that EndoWorks supports.

Windows Updates	Description	Tested	Pass/Fail	Date
KB890830	Windows Malicious Software Removal Tool - July 2015	Yes	Pass	7/17/2015
KB2934207	Windows XP End of Support Notification	Yes	Pass	7/17/2015
KB975558	Security Update for Windows XP	Yes	Pass	7/17/2015
KB954155	Security Update for Windows Media Format Runtime 9.9.5811 for Windows XP SP3	Yes	Pass	7/17/2015
KB953300	Microsoft .NET Framework 2.0 Service Pack 1 Security Update for Windows 2000, Windows Server 2003, and Windows XP	Yes	Pass	7/17/2015
KB905474	Windows Genuine Advantage Notification	Yes	Pass	7/17/2015
KB941569	Security Update for Windows XP with Windows Media Format Runtime 9.5 and 11	Yes	Pass	7/17/2015
KB892130	Windows Genuine Advantage Validation Tool	Yes	Pass	7/17/2015
KB951847	Microsoft .NET Framework 3.5 Service Pack 1 and .Net Framework 3.5 Family Update for .Net versions 2.0 through 3.5 x86	Yes	Pass	7/17/2015
KB2964358	Security update for Internet Explorer 8 for Windows XP	Yes	Pass	5/2/2014
KB2936068	Cumulative Security Update for Internet Explorer 8 for Windows XP	Yes	Pass	4/24/2014
KB890830	Windows Malicious Software Removal Tool - April 2014	Yes	Pass	4/24/2014
KB2922229	Security Update for Windows XP	Yes	Pass	4/24/2014
KB2925418	Cumulative Security Update for Internet Explorer 8 for Windows XP	Yes	Pass	3/24/2014
KB2929961	Security Update for Windows XP	Yes	Pass	3/24/2014
KB2930275	Security Update for Windows XP	Yes	Pass	3/24/2014
KB2934207	Windows XP End of Support Notification	Yes	Pass	3/24/2014
KB890830	Windows Malicious Software Removal Tool - March 2014	Yes	Pass	3/24/2014
KB2836941	Update for Microsoft .NET Framework 2.0 SP2 on Windows Server 2003 and Windows XP x86	Yes	Pass	2/12/2014
KB2898856	Security update for Microsoft.NET Framework 2.0 SP2 on Window Server 2003 and Windows XP x86	Yes	Pass	2/12/2014
KB2901111	Security update for Microsoft.NET Framework 2.0 SP2 on Window Server 2003 and Windows XP x86	Yes	Pass	2/12/2014
KB2909210	Security Update for Windows XP	Yes	Pass	2/12/2014
KB2909921	Cumulative Security Update for Internet Explorer 8 for Windows XP	Yes	Pass	2/12/2014
KB2916036	Security Update for Windows XP	Yes	Pass	2/12/2014
KB890830	Windows Malicious Software Removal Tool - February 2014	Yes	Pass	2/12/2014
KB2909212	Security Update for Windows XP	Yes	Pass	2/12/2014
KB2510531	Security Update for Windows XP	Yes	Pass	1/23/2014
KB2898785	Cumulative Security Update for Internet Explorer 8 for Windows XP	Yes	Pass	1/23/2014
KB2914368	Security Update for Windows XP	Yes	Pass	1/23/2014
KB890830	Windows Malicious Software Removal Tool - January 2014	Yes	Pass	1/23/2014
KB2892075	Security Update for Windows XP	Yes	Pass	12/18/2013
KB2893294	Security Update for Windows XP	Yes	Pass	12/18/2013
KB2893984	Security Update for Windows XP	Yes	Pass	12/18/2013
KB2898715	Security Update for Windows XP	Yes	Pass	12/18/2013
KB2898785	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	12/18/2013
KB2904266	Update for Windows XP	Yes	Pass	12/18/2013
KB2917500	Security Update for Windows XP Windows Server 2003	Yes	Pass	12/18/2013
KB890830	Windows Malicious Software Removal Tool - December 2013	Yes	Pass	12/18/2013
KB2862152	Security Update for Windows XP	Yes	Pass	11/21/2013
KB2868626	Security Update for Windows XP	Yes	Pass	11/21/2013
KB2876331	Security Update for Windows XP	Yes	Pass	11/21/2013
KB2888505	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	11/21/2013
KB2900986	Cumulative Security Update for ActiveX Killbits for Windows XP	Yes	Pass	11/21/2013
KB890830	Windows Malicious Software Removal Tool - November 2013	Yes	Pass	11/21/2013
KB2847311	Security Update for Windows XP	Yes	Pass	10/23/2013
KB2861189	Security Update for Microsoft .NET Framework 3.0 SP2 on Windows Server 2003 and Windows XP x86	Yes	Pass	10/23/2013
KB2861697	Security Update for Microsoft .NET Framework 3.5 SP1 on Windows XP, Server 2003, Vista, Server 2008 x86	Yes	Pass	10/23/2013
KB2862330	Security Update for Windows XP	Yes	Pass	10/23/2013
KB2862335	Security Update for Windows XP	Yes	Pass	10/23/2013
KB2863239	Security Update for Microsoft .NET Framework 2.0 SP2 on Windows Server 2003 and Windows XP x86	Yes	Pass	10/23/2013
KB2879017	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	10/23/2013
KB2883150	Security Update for Windows XP	Yes	Pass	10/23/2013
KB2884256	Security Update for Windows XP	Yes	Pass	10/23/2013
KB890830	Windows Malicious Software Removal Tool - October 2013	Yes	Pass	10/23/2013
KB2834903	Security Update for Windows Media Format Runtime 9.5 for Windows XP	Yes	Pass	9/23/2013
KB2864063	Security Update for Windows XP	Yes	Pass	9/23/2013
KB2870699	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	9/23/2013
KB2876217	Security Update for Windows XP	Yes	Pass	9/23/2013
KB2876315	Security Update for Windows XP	Yes	Pass	9/23/2013
KB890830	Windows Malicious Software Removal Tool - September 2013	Yes	Pass	9/23/2013
KB2834903	Security Update for Windows Media Format Runtime 9.5 for Windows XP	Yes	Pass	8/21/2013
KB2844285	Security Update for Microsoft .NET Framework 2.0 SP2 on Windows Server 2003 and Windows XP x86	Yes	Pass	8/21/2013
KB2849470	Security Update for Windows XP	Yes	Pass	8/21/2013
KB2850869	Security Update for Windows XP	Yes	Pass	8/21/2013
KB2859537	Security Update for Windows XP	Yes	Pass	8/21/2013
KB2862772	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	8/21/2013
KB2863058	Update for windows XP	Yes	Pass	8/21/2013
KB890830	Windows Malicious Software Removal Tool - August 2013	Yes	Pass	8/21/2013
KB2832411	Security Update for Microsoft .NET Framework 3.0 SP2 on Windows Server 2003 and Windows XP x86	Yes	Pass	7/16/2013

KB2833940	Security Update for Microsoft .NET Framework 2.0 SP2 on Windows Server 2003 and Windows XP x86	Yes	Pass	7/16/2013
KB2833941	Security Update for Microsoft .NET Framework 1.1 SP1 on Windows XP, Windows Vista and Windows Server 2008 x86	Yes	Pass	7/16/2013
KB2834886	Security Update for Windows XP	Yes	Pass	7/16/2013
KB2834903	Security Update for Windows Media Format Runtime 9.5 for Windows XP	Yes	Pass	7/16/2013
KB2840629	Security Update for Microsoft .NET Framework 3.5 SP1 on Windows XP, Server 2003, Vista, Server 2008 x86	Yes	Pass	7/16/2013
KB2844285	Security Update for Microsoft .NET Framework 2.0 SP2 on Windows Server 2003 and Windows XP x86	Yes	Pass	7/16/2013
KB2845187	Security Update for Windows XP	Yes	Pass	7/16/2013
KB2846071	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	7/16/2013
KB2850851	Security Update for Windows XP	Yes	Pass	7/16/2013
KB890830	Windows Malicious Software Removal Tool - July 2013	Yes	Pass	7/16/2013
KB2838727	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	6/12/2013
KB2839229	Security Update for Windows XP	Yes	Pass	6/12/2013
KB890830	Windows Malicious Software Removal Tool - June 2013	Yes	Pass	6/12/2013
KB2804577	Security Update for Microsoft .NET Framework 2.0 SP2 on Windows Server 2003 and Windows XP x86	Yes	Pass	5/17/2013
KB2820197	Cumulative Security Update for ActiveX Killbits for Windows XP	Yes	Pass	5/17/2013
KB2829361	Security Update for Windows XP	Yes	Pass	5/17/2013
KB2829530	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	5/17/2013
KB890830	Windows Malicious Software Removal Tool - May 2013	Yes	Pass	5/17/2013
KB2808735	Security Update for Windows XP	Yes	Pass	4/23/2013
KB2813170	Security Update for Windows XP	Yes	Pass	4/23/2013
KB2813345	Security Update for Windows XP	Yes	Pass	4/23/2013
KB2817183	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	4/23/2013
KB2820917	Security Update for Windows XP	Yes	Pass	4/23/2013
KB890830	Windows Malicious Software Removal Tool - April 2013	Yes	Pass	4/23/2013
KB2807986	Security Update for Windows XP	Yes	Pass	3/23/2013
KB2809289	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	3/23/2013
KB890830	Windows Malicious Software Removal Tool - March 2013	Yes	Pass	3/23/2013
KB2778344	Security Update for Windows XP	Yes	Pass	3/6/2013
KB2780091	Security Update for Windows XP	Yes	Pass	3/6/2013
KB2789643	Security Update for Microsoft .NET Framework 2.0 SP2 on Windows Server 2003 and Windows XP x86	Yes	Pass	3/6/2013
KB2792100	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	3/6/2013
KB2797052	Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	3/6/2013
KB2799494	Security Update for Windows XP	Yes	Pass	3/6/2013
KB2802968	Security Update for Windows XP	Yes	Pass	3/6/2013
KB890830	Windows Malicious Software Removal Tool - February 2013	Yes	Pass	3/6/2013
KB2736416	Security Update for Microsoft .NET Framework 3.5 SP1 on Windows XP, Server 2003, Vista, Server 2008 x86	Yes	Pass	1/17/2013
KB2742596	Security Update for Microsoft .NET Framework 2.0 SP2 on Windows Server 2003 and Windows XP x86	Yes	Pass	1/17/2013
KB2742597	Security Update for Microsoft .NET Framework 1.1 SP1 on Windows XP, Windows Vista, and Windows Server 2008 x86	Yes	Pass	1/17/2013
KB2753842	Security Update for Windows XP	Yes	Pass	1/17/2013
KB2756918	Security Update for Microsoft .NET Framework 3.0 SP2 on Windows Server 2003 and Windows XP x86	Yes	Pass	1/17/2013
KB2757638	Security Update for Windows XP	Yes	Pass	1/17/2013
KB2798897	Update for Windows XP and Windows Server 2003	Yes	Pass	1/17/2013
KB2799329	Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	1/17/2013
KB890830	Windows Malicious Software Removal Tool - January 2013	Yes	Pass	1/17/2013
KB2753842	Security Update for Windows XP	Yes	Pass	12/28/2012
KB2758857	Security Update for Windows XP	Yes	Pass	12/28/2012
KB2761465	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	12/28/2012
KB2770660	Security Update for Windows XP	Yes	Pass	12/28/2012
KB2779030	Security Update for Windows XP	Yes	Pass	12/28/2012
KB2779562	Update for Windows XP	Yes	Pass	12/28/2012
KB890830	Windows Malicious Software Removal Tool - December 2012	Yes	Pass	12/28/2012
KB2698023	Security Update for Microsoft .NET Framework 1.1 SP1 on Windows XP, Windows Vista, and Windows Server 2008 x86	Yes	Pass	11/19/2012
KB2727528	Security Update for Windows XP	Yes	Pass	11/19/2012
KB2729450	Security Update for Microsoft .NET Framework 2.0 SP2 on Windows Server 2003 and Windows XP x86	Yes	Pass	11/19/2012
KB2761226	Security Update for Windows XP	Yes	Pass	11/19/2012
KB890830	Windows Malicious Software Removal Tool - November 2012	Yes	Pass	11/19/2012
KB2661254	Update for Windows XP	Yes	Pass	10/17/2012
KB2724197	Security Update for Windows XP	Yes	Pass	10/17/2012
KB2744842	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	10/17/2012
KB2749655	Update for Windows XP	Yes	Pass	10/17/2012
KB2756822	Update for Windows XP	Yes	Pass	10/17/2012
KB890830	Windows Malicious Software Removal Tool - October 2012	Yes	Pass	10/17/2012
KB2736233	Update Rollup for ActiveX Killbits for Windows XP	Yes	Pass	9/18/2012
KB890830	Windows Malicious Software Removal Tool - September 2012	Yes	Pass	9/18/2012
KB2705219	Security Update for Windows XP	Yes	Pass	8/22/2012
KB2712808	Security Update for Windows XP	Yes	Pass	8/22/2012
KB2722913	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	8/22/2012
KB2723135	Security Update for Windows XP	Yes	Pass	8/22/2012
KB2731847	Security Update for Windows XP	Yes	Pass	8/22/2012
KB890830	Windows Malicious Software Removal Tool - August 2012	Yes	Pass	8/22/2012
KB2655992	Security Update for Windows XP	Yes	Pass	7/18/2012
KB2691442	Security Update for Windows XP	Yes	Pass	7/18/2012
KB2698365	Security Update for Windows XP	Yes	Pass	7/18/2012
KB2718523	Security Update for Windows XP	Yes	Pass	7/18/2012
KB2719985	Security Update for Windows XP	Yes	Pass	7/18/2012
KB2728973	Update for Windows XP and Windows Server 2003	Yes	Pass	7/18/2012
KB890830	Windows Malicious Software Removal Tool - July 2012	Yes	Pass	7/18/2012
KB2656369	Security Update for Microsoft .NET Framework 2.0 SP2 on Windows Server 2003 and Windows XP x86	Yes	Pass	6/21/2012
KB2656370	Security Update for Microsoft .NET Framework 1.1 SP1 on Windows XP, Windows Vista, and Windows Server 2008 x86	Yes	Pass	6/21/2012
KB2685939	Security Update for Windows XP	Yes	Pass	6/21/2012
KB2688828	Security Update for Microsoft .NET Framework 2.0 SP2 on Windows Server 2003 and Windows XP x86	Yes	Pass	6/21/2012
KB2699988	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	6/21/2012
KB2707511	Security Update for Windows XP	Yes	Pass	6/21/2012
KB2709162	Security Update for Windows XP	Yes	Pass	6/21/2012
KB2718704	Update for Windows XP	Yes	Pass	6/21/2012
KB890830	Windows Malicious Software Removal Tool - June 2012	Yes	Pass	6/21/2012
KB2604092	Security Update for Microsoft .NET Framework 2.0 SP2 on Windows Server 2003 and Windows XP x86	Yes	Pass	5/10/2012
KB2604110	Security Update for Microsoft .NET Framework 3.0 SP2 on Windows Server 2003 and Windows XP x86	Yes	Pass	5/10/2012
KB2604111	Security Update for Microsoft .NET Framework 3.5 SP1 on Windows XP, Server 2003, Vista, Server 2008 x86	Yes	Pass	5/10/2012

KB2656407	Security Update for Microsoft .NET Framework 3.0 SP2 on Windows Server 2003 and Windows XP x86	Yes	Pass	5/10/2012
KB2659262	Security Update for Windows XP	Yes	Pass	5/10/2012
KB2676562	Security Update for Windows XP	Yes	Pass	5/10/2012
KB2686509	Security Update for Windows XP	Yes	Pass	5/10/2012
KB2695962	Update Rollup for ActiveX Killbits for Windows XP	Yes	Pass	5/10/2012
KB890830	Windows Malicious Software Removal Tool - May 2012	Yes	Pass	5/10/2012
KB2653956	Security Update for Windows XP	Yes	Pass	4/16/2012
KB2660649	Security Update for Windows XP	Yes	Pass	4/16/2012
KB2656369	Security Update for Microsoft .NET Framework 2.0 SP2 on Windows Server 2003 and Windows XP x86	Yes	Pass	4/16/2012
KB2656370	Security Update for Microsoft .NET Framework 1.1 SP1 on Windows XP, Windows Vista and Windows Server 2008 x86	Yes	Pass	4/16/2012
KB2675157	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	4/16/2012
KB890830	Windows Malicious Software Removal Tool - April 2012	Yes	Pass	4/16/2012
KB2621440	Security Update for Windows XP	Yes	Pass	3/16/2012
KB2641653	Security Update for Windows XP	Yes	Pass	3/16/2012
KB2647518	Update Rollup for Active X Killbits for Windows XP	Yes	Pass	3/16/2012
KB890830	Windows Malicious Software Removal Tool - March 2012	Yes	Pass	3/16/2012
KB2633880	Security Update for Microsoft .NET Framework 2.0 SP2 on Windows Server 2003 and Windows XP x86	Yes	Pass	2/17/2012
KB2647516	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	2/17/2012
KB2660465	Security Update for Windows XP	Yes	Pass	2/17/2012
KB2661637	Security Update for Windows XP	Yes	Pass	2/17/2012
KB890830	Windows Malicious Software Removal Tool - February 2012	Yes	Pass	2/17/2012
KB2584146	Security Update for Windows XP	Yes	Pass	1/18/2012
KB2585542	Security Update for Windows XP	Yes	Pass	1/18/2012
KB2598479	Security Update for Windows XP	Yes	Pass	1/18/2012
KB2603381	Security Update for Windows XP	Yes	Pass	1/18/2012
KB2631813	Security Update for Windows XP	Yes	Pass	1/18/2012
KB2646524	Security Update for Windows XP	Yes	Pass	1/18/2012
KB2656352	Security Update for Microsoft .NET Framework 2.0 SP2 on Windows Server 2003 and Windows XP x86	Yes	Pass	1/18/2012
KB2656353	Security Update for Microsoft .NET Framework 1.1 SP1 on Windows XP, Windows Vista and Windows Server 2008 x86	Yes	Pass	1/18/2012
KB2632503	Update for Windows XP	Yes	Pass	1/18/2012
KB2657424	Security Update for Microsoft .NET Framework 3.5 SP1 on Windows XP, Server 2003, Vista, Server 2008 x86	Yes	Pass	1/18/2012
KB890830	Windows Malicious Software Removal Tool - January 2012	Yes	Pass	1/18/2012
KB2618444	Cumulative Security Update for Internet Explorer for Windows XP	Yes	Pass	12/26/2011
KB2618451	Cumulative Security Update for ActiveX Killbits for Windows XP	Yes	Pass	12/26/2011
KB2619339	Security Update for Windows XP	Yes	Pass	12/26/2011
KB2620712	Security Update for Windows XP	Yes	Pass	12/26/2011
KB2624667	Security Update for Windows XP	Yes	Pass	12/26/2011
KB2633171	Security Update for Windows XP	Yes	Pass	12/26/2011
KB2633952	Update for Windows XP	Yes	Pass	12/26/2011
KB2639417	Security Update for Windows XP	Yes	Pass	12/26/2011
KB890830	Windows Malicious software Removal Tool-December 2011	Yes	Pass	12/26/2011
KB2544893	Security Update for Windows XP	Yes	Pass	11/15/2011
KB2641690	Update for Windows XP	Yes	Pass	11/15/2011
KB890830	Windows Malicious Software Removal Tool - November 2011	Yes	Pass	11/15/2011
KB2518864	Security update for .Net Framework 2.0 SP2 and 3.5 SP1 on Windows Server 2003 And Windows XP x86	Yes	Pass	10/29/2011
KB2564958	Security Update for Windows XP	Yes	Pass	10/29/2011
KB2567053	Security Update for Windows XP	Yes	Pass	10/29/2011
KB2572067	Security Update for Microsoft .Net Framework 1.1 SP1 on Windows XP,Vista 2003 and Windows Server 2008 x86	Yes	Pass	10/29/2011
KB2572073	Security update for .Net Framework 2.0 SP2 and 3.5 SP1 on Windows Server 2003 and Windows XP x86	Yes	Pass	10/29/2011
KB2586448	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	10/29/2011
KB2592799	Security Update for Windows XP	Yes	Pass	10/29/2011
KB890830	Windows Malicious Software Removal Tool -October 2011	Yes	Pass	10/29/2011
KB961118	Update for Windows XP	Yes	Pass	10/29/2011
KB963707	Update to .NET Framework 3.5 Service Pack 1 for the .NET Framework assistant 1.0 x86	Yes	Pass	10/29/2011
KB2570791	Update for Windows XP	Yes	Pass	9/22/2011
KB2570947	Security Update for Windows XP	Yes	Pass	9/22/2011
KB2616676	Update for Windows XP	Yes	Pass	9/22/2011
KB890830	Windows Malicious Software Removal Tool - September 2011	Yes	Pass	9/22/2011
KB2536276	Security Update for Windows XP	Yes	Pass	8/18/2011
KB2539631	Security Update for Microsoft .Net Framework 2.0 SP2 on Windows Server 2003 and Windows XP x86	Yes	Pass	8/18/2011
KB2559049	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	8/18/2011
KB2562937	Update Rollup for Active X Killbits for Windows XP	Yes	Pass	8/18/2011
KB2566454	Security Update for Windows XP	Yes	Pass	8/18/2011
KB2567680	Security Update for Windows XP	Yes	Pass	8/18/2011
KB2570222	Security Update for Windows XP	Yes	Pass	8/18/2011
KB890830	Windows Malicious Software Removal Tool-August 2011	Yes	Pass	8/18/2011
KB2507938	Security Update for Windows XP	Yes	Pass	7/18/2011
KB2541763	Update for Windows XP	Yes	Pass	7/18/2011
KB2555917	Security Update for Windows XP	Yes	Pass	7/18/2011
KB890830	Windows Malicious Software Removal Tool-July 2011	Yes	Pass	7/18/2011
KB2476490	Security Update for Windows XP	Yes	Pass	6/21/2011
KB2478658	and Windows XP x86	Yes	Pass	6/21/2011
KB2503665	Security Update for Windows XP	Yes	Pass	6/21/2011
KB2518864	Security update for .Net Framework 2.0 SP2 and 3.5 SP1 on Windows Server 2003 and Windows XP x86	Yes	Pass	6/21/2011
KB2530548	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	6/21/2011
KB2535512	Security Update for Windows XP	Yes	Pass	6/21/2011
KB2536276	Security Update for Windows XP	Yes	Pass	6/21/2011
KB2544521	Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	6/21/2011
KB2483614	Security Update for Windows XP	Yes	Pass	6/21/2011
KB2544893	Security Update for Windows XP	Yes	Pass	6/21/2011
KB890830	Windows Malicious Software Removal Tool-June 2011	Yes	Pass	6/21/2011
KB890830	Windows Malicious Software Removal Tool-May 2011	Yes	Pass	5/20/2011
KB2412687	Security Update for Windows XP	Yes	Pass	4/14/2011
KB2485663	Security Update for Windows XP	Yes	Pass	4/14/2011
KB2491683	Security Update for Windows XP	Yes	Pass	4/14/2011
KB2497640	Cumulative Security Update for Internet Explorer 6 fr Windows XP	Yes	Pass	4/14/2011
KB2503658	Security Update for Windows XP	Yes	Pass	4/14/2011
KB2506212	Security Update for Windows XP	Yes	Pass	4/14/2011
KB2506223	Security Update for Windows XP	Yes	Pass	4/14/2011

KB2507618	Security Update for Windows XP	Yes	Pass	4/14/2011
KB2508272	Cumulative Security Update for ActiveX Killbits for Windows XP	Yes	Pass	4/14/2011
KB2508429	Security Update for Windows XP	Yes	Pass	4/14/2011
KB2509553	Security Update for Windows XP	Yes	Pass	4/14/2011
KB2510581	Security Update for Windows XP	Yes	Pass	4/14/2011
KB2511455	Security Update for Windows XP	Yes	Pass	4/14/2011
KB2524375	Update for Windows XP	Yes	Pass	4/14/2011
KB890830	Windows Malicious Software Removal Tool-April 2011	Yes	Pass	4/14/2011
KB2479943	Security Update for Windows XP	Yes	Pass	3/15/2011
KB2481109	Security Update for Windows XP	Yes	Pass	3/15/2011
KB890830	Windows Malicious Software Removal Tool - March 2011	Yes	Pass	3/15/2011
KB971029	Update for Windows XP	Yes	Pass	3/15/2011
KB2393802	Security Update for Windows XP	Yes	Pass	2/22/2011
KB2476687	Security Update for Windows XP	Yes	Pass	2/22/2011
KB2478960	Security Update for Windows XP	Yes	Pass	2/22/2011
KB2478971	Security Update for Windows XP	Yes	Pass	2/22/2011
KB2479628	Security Update for Windows XP	Yes	Pass	2/22/2011
KB2482017	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	2/22/2011
KB2483185	Security Update for Windows XP	Yes	Pass	2/22/2011
KB2485376	Security Update for Windows XP	Yes	Pass	2/22/2011
KB890830	Windows Malicious Software Removal Tool - February 2011	Yes	Pass	2/22/2011
KB2419632	Security Update for Windows XP	Yes	Pass	1/24/2011
KB890830	Windows Malicious Software Removal Tool - January 2011	Yes	Pass	1/24/2011
KB2296199	Security Update for Windows XP	Yes	Pass	12/27/2010
KB2416400	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	12/27/2010
KB2423089	Security Update for Windows XP	Yes	Pass	12/27/2010
KB2436673	Security Update for Windows XP	Yes	Pass	12/27/2010
KB2440591	Security Update for Windows XP	Yes	Pass	12/27/2010
KB2443105	Security Update for Windows XP	Yes	Pass	12/27/2010
KB2443685	Update for Windows XP	Yes	Pass	12/27/2010
KB2467659	Update for Internet Explorer for Windows XP	Yes	Pass	12/27/2010
KB890830	Windows Malicious Software Removal Tool - December 2010	Yes	Pass	12/27/2010
KB2345009	Security Update for Microsoft Word Viewer	Yes	Pass	11/29/2010
KB890830	Windows Malicious Software Removal Tool - November 2010	Yes	Pass	11/29/2010
KB2158563	Update for Windows XP	Yes	Pass	10/27/2010
KB2279986	Security Update for Windows XP	Yes	Pass	10/27/2010
KB2296011	Security Update for Windows XP	Yes	Pass	10/27/2010
KB2345886	Update for Windows XP	Yes	Pass	10/27/2010
KB2360131	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	10/27/2010
KB2360937	Security Update for Windows XP	Yes	Pass	10/27/2010
KB2378111	Security Update for Windows XP	Yes	Pass	10/27/2010
KB2387149	Security Update for Windows XP	Yes	Pass	10/27/2010
KB2416447	Security Update for Microsoft .Net Framework 1.1 SP1 on Windows XP, Windows Vista, and Windows Server 2008 x86	Yes	Pass	10/27/2010
KB2416473	Security Update for Microsoft .Net Framework 3.5 SP1 Windows XP, Windows Server 2003, Windows Vista, Windows Server 2008 x86	Yes	Pass	10/27/2010
KB2418241	Security Update for Microsoft .Net Framework 2.0 SP2 and 3.5 SP1 on Windows Server 2003, Windows XP x86	Yes	Pass	10/27/2010
KB890830	Windows Malicious Software Removal Tool - October 2010	Yes	Pass	10/27/2010
KB979687	Security Update for Windows XP	Yes	Pass	10/27/2010
KB981957	Security Update for Windows XP	Yes	Pass	10/27/2010
KB982132	Security Update for Windows XP	Yes	Pass	10/27/2010
KB2121546	Security Update for Windows XP	Yes	Pass	9/30/2010
KB2141007	Update for Windows XP	Yes	Pass	9/30/2010
KB2259922	Security Update for Windows XP	Yes	Pass	9/30/2010
KB2347290	Security Update for Windows XP	Yes	Pass	9/30/2010
KB890830	Windows Malicious Software Removal Tool - September 2010	Yes	Pass	9/30/2010
KB975558	Security Update for Windows XP	Yes	Pass	9/30/2010
KB981322	Security Update for Windows XP	Yes	Pass	9/30/2010
KB982802	Security Update for Windows XP	Yes	Pass	9/30/2010
KB2079403	Security Update for Windows XP	Yes	Pass	8/12/2010
KB2115168	Security Update for Windows XP	Yes	Pass	8/12/2010
KB2160329	Security Update for Windows XP	Yes	Pass	8/12/2010
KB2183461	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	8/12/2010
KB2251437	Security Update for Microsoft Word Viewer	Yes	Pass	8/12/2010
KB2286198	Security Update for Windows XP	Yes	Pass	8/12/2010
KB890830	Windows Malicious Software Removal Tool - August 2010	Yes	Pass	8/12/2010
KB980436	Security Update for Windows XP	Yes	Pass	8/12/2010
KB981852	Security Update for Windows XP	Yes	Pass	8/12/2010
KB981997	Security Update for Windows XP	Yes	Pass	8/12/2010
KB982214	Security Update for Windows XP	Yes	Pass	8/12/2010
KB982665	Security Update for Windows XP	Yes	Pass	8/12/2010
KB983583	Security update for .Net Framework 2.0 SP2 and 3.5 SP1 on Windows Server 2003 and Windows XP x86	Yes	Pass	8/12/2010
KB2229593	Security Update for Windows XP	Yes	Pass	7/27/2010
KB890830	Windows Malicious Software Removal Tool - July 2010	Yes	Pass	7/27/2010
KB890830	Windows Malicious Software Removal Tool - June 2010	Yes	Pass	6/24/2010
KB975562	Security Update for Windows XP	Yes	Pass	6/24/2010
KB978695	Security Update for Windows Media Format Runtime 9.9.5&11 for Windows XP SP3	Yes	Pass	6/24/2010
KB979482	Security Update for Windows XP	Yes	Pass	6/24/2010
KB979559	Security Update for Windows XP	Yes	Pass	6/24/2010
KB979906	Microsoft .Net Framework 1.1 SP1 Security Update for Windows 2000, Windows XP	Yes	Pass	6/24/2010
KB979909	Microsoft .Net Framework 3.5 SP1 and .Net Framework 2.0 SP2 Security Update for Windows 2000, Windows Server 2003, Windows XP x86	Yes	Pass	6/24/2010
KB980195	Cumulative Security Update for ActiveX Killbits for Windows XP	Yes	Pass	6/24/2010
KB980218	Security Update for Windows XP	Yes	Pass	6/24/2010
KB982168	Microsoft .Net Framework 3.5 SP1 Update for Windows Server 2003, Windows XP x86	Yes	Pass	6/24/2010
KB982381	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	6/24/2010
KB982524	Microsoft .Net Framework 3.5 SP1 and .Net Framework 2.0 SP2 Update for Windows Server 2003, and Windows XP x86	Yes	Pass	6/24/2010
KB890830	Windows Malicious Software Removal Tool - May 2010	Yes	Pass	5/25/2010
KB978542	Security Update for Windows XP	Yes	Pass	5/25/2010
KB890830	Windows Malicious Software Removal Tool - April 2010	Yes	Pass	4/27/2010
KB977816	Security Update for Windows XP	Yes	Pass	4/27/2010
KB978338	Security Update for Windows XP	Yes	Pass	4/27/2010

KB978601	Security Update for Windows XP	Yes	Pass	4/27/2010
KB979309	Security Update for Windows XP	Yes	Pass	4/27/2010
KB979402	Security Update for Windows Media Player	Yes	Pass	4/27/2010
KB979683	Security Update for Windows XP	Yes	Pass	4/27/2010
KB980182	Update for Windows XP	Yes	Pass	4/27/2010
KB980232	Security Update for Windows XP	Yes	Pass	4/27/2010
KB981349	Security Update for Windows XP	Yes	Pass	4/27/2010
KB890830	Windows Malicious Software Removal Tool - March 2010	Yes	Pass	3/16/2010
KB975561	Security Update for Windows XP	Yes	Pass	3/16/2010
KB979306	Update for Windows XP	Yes	Pass	3/16/2010
KB890830	Windows Malicious software Removal Tool - February 2010	Yes	Pass	2/22/2010
KB971468	Security Update for Windows XP	Yes	Pass	2/22/2010
KB975560	Security Update for Windows XP	Yes	Pass	2/22/2010
KB975713	Security Update for Windows XP	Yes	Pass	2/22/2010
KB977165	Security Update for Windows XP	Yes	Pass	2/22/2010
KB977914	Security Update for Windows XP	Yes	Pass	2/22/2010
KB978037	Security Update for Windows XP	Yes	Pass	2/22/2010
KB978207	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	2/22/2010
KB978251	Security Update for Windows XP	Yes	Pass	2/22/2010
KB978262	Cumulative Security Update for ActiveX Killbits for Windows XP	Yes	Pass	2/22/2010
KB978706	Security Update for Windows XP	Yes	Pass	2/22/2010
KB890830	Windows Malicious Software Removal Tool - January 2010	Yes	Pass	1/21/2010
KB972270	Security Update for Windows Server 2003	Yes	Pass	1/21/2010
KB978557	Update for Microsoft Office Excel Viewer 2003	Yes	Pass	1/21/2010
KB978558	Update for Microsoft Word Viewer 2003	Yes	Pass	1/21/2010
KB890830	Windows Malicious Software Removal Tool - December 2009	Yes	Pass	12/22/2009
KB955759	Update for Windows XP	Yes	Pass	12/22/2009
KB970430	Update for Windows XP	Yes	Pass	12/22/2009
KB971737	Update for Windows XP	Yes	Pass	12/22/2009
KB973484	Security Update for Microsoft Excel Viewer 2003	Yes	Pass	12/22/2009
KB973866	Security Update for Microsoft Word Viewer 2003	Yes	Pass	12/22/2009
KB973904	Security Update for Windows XP	Yes	Pass	12/22/2009
KB974318	Security Update for Windows XP	Yes	Pass	12/22/2009
KB974392	Security Update for Windows XP	Yes	Pass	12/22/2009
KB976325	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	12/22/2009
KB890830	Windows Malicious Software Removal Tool - November 2009	Yes	Pass	12/7/2009
KB949810	Office Genuine Advantage Notifications	Yes	Pass	12/7/2009
KB969947	Security Update for Windows XP	Yes	Pass	12/7/2009
KB973687	Security Update for Windows - November 2009	Yes	Pass	12/7/2009
KB973688	Update for Microsoft XML core Services 4.0 Service Pack 2	Yes	Pass	12/7/2009
KB976098	Update for Windows XP	Yes	Pass	12/7/2009
KB976749	Update for Internet Explorer 6 for Windows XP	Yes	Pass	12/7/2009
KB890830	Windows Malicious Software Removal Tool - October 2009	Yes	Pass	10/24/2009
KB953155	Security Update for Windows XP	Yes	Pass	10/24/2009
KB958869	Security Update for Windows XP	Yes	Pass	10/24/2009
KB969059	Security Update for Windows XP	Yes	Pass	10/24/2009
KB971486	Security Update for Windows XP	Yes	Pass	10/24/2009
KB974112	Security Update for Windows XP	Yes	Pass	10/24/2009
KB974455	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	10/24/2009
KB974571	Security Update for Windows XP	Yes	Pass	10/24/2009
KB975025	Security Update for Windows XP	Yes	Pass	10/24/2009
KB890830	Windows Malicious Software Removal Tool - September 2009	Yes	Pass	9/14/2009
KB956844	Security Update for Windows XP	Yes	Pass	9/14/2009
KB968816	Security Update for Windows Media Format Runtime 9.9.5&11 for Windows XP SP3	Yes	Pass	9/14/2009
KB970653	Update for Windows XP	Yes	Pass	9/14/2009
KB971961	Security Update for Jscript 5.7 for Windows XP	Yes	Pass	9/14/2009
KB890830	Windows Malicious Software Removal Tool - August 2009	Yes	Pass	8/17/2009
KB905474	Windows Genuine Advantage Notification	Yes	Pass	8/17/2009
KB923561	Security Update for Windows XP	Yes	Pass	8/17/2009
KB951847	Microsoft .Net Framework 3.5 service pack 1 and .Net Framework 3.5 Family Update for .Net versions 2.0 through 3.5 X86	Yes	Pass	8/17/2009
KB952004	Security Update for Windows XP	Yes	Pass	8/17/2009
KB952069	Security Update for Windows XP Service Pack 3	Yes	Pass	8/17/2009
KB954430	Security Update for Windows XP	Yes	Pass	8/17/2009
KB954459	Security Update for Windows XP	Yes	Pass	8/17/2009
KB954600	Security Update for Windows XP	Yes	Pass	8/17/2009
KB955069	Security Update for Windows XP	Yes	Pass	8/17/2009
KB955839	Update for Windows XP	Yes	Pass	8/17/2009
KB956572	Security Update for Windows XP	Yes	Pass	8/17/2009
KB956744	Security Update for Windows XP	Yes	Pass	8/17/2009
KB956802	Security Update for Windows XP	Yes	Pass	8/17/2009
KB957097	Security Update for Windows XP	Yes	Pass	8/17/2009
KB958644	Security Update for Windows XP (Critical Update)	Yes	Pass	8/17/2009
KB958687	Security Update for Windows XP	Yes	Pass	8/17/2009
KB959426	Security Update for Windows XP	Yes	Pass	8/17/2009
KB960225	Security Update for Windows XP	Yes	Pass	8/17/2009
KB960803	Security Update for Windows XP	Yes	Pass	8/17/2009
KB960859	Security Update for Windows XP	Yes	Pass	8/17/2009
KB961371	Security Update for Windows XP	Yes	Pass	8/17/2009
KB961501	Security Update for Windows XP	Yes	Pass	8/17/2009
KB967715	Update for Windows XP	Yes	Pass	8/17/2009
KB968389	Update for Windows XP	Yes	Pass	8/17/2009
KB968537	Security Update for Windows XP	Yes	Pass	8/17/2009
KB970238	Security Update for Windows XP	Yes	Pass	8/17/2009
KB971557	Security Update for Windows XP	Yes	Pass	8/17/2009
KB971633	Security Update for Windows XP	Yes	Pass	8/17/2009
KB971657	Security Update for Windows XP	Yes	Pass	8/17/2009
KB972260	Cumulative Security Update for Internet Explorer 6 for Windows XP	Yes	Pass	8/17/2009
KB973346	Cumulative Security Update for ActiveX Killbits for Windows XP	Yes	Pass	8/17/2009
KB973354	Security Update for Windows XP	Yes	Pass	8/17/2009

KB973507	Security Update for Windows XP	Yes	Pass	8/17/2009
KB973540	Security Update for Windows XP Service Pack 3	Yes	Pass	8/17/2009
KB973815	Security Update for Windows XP	Yes	Pass	8/17/2009
KB973869	Security Update for Windows XP	Yes	Pass	8/17/2009